

**FOREIGN BODY PROTECTION CAP -SINGLE**

END-WS0420000

**PRODUCT DESCRIPTION**

**Designed to fit on the distal end of your endoscope to form a protective "shield" (or hood) around your biopsy sample or foreign body retrieval while retracting the endoscope from the patient**

Will fit on most gastroscopes and colonoscopes - a thick rubber band supplied aids in fitting on smaller endoscopes or when larger diameter endoscopes are used since they require the use of a lubricant to install without damaging the distal tip. Hoods come individually packed in a box of ten (10) hoods with instructions for use.

*Devices are made of natural rubber latex. Those that are allergic to latex should not handle or ingest these hoods and they should not be used in patients that are allergic to latex.*

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	END-WS0420000
<b>Procedure Type</b>	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights