

## BURTONS NON VACUUM VETCLAVE

2545-V S3 A



### PRODUCT DESCRIPTION

## The Vetclave 2545 non-vacuum series 3 model has been designed specifically to meet the needs of your practice.

You have the ability to increase the sterilisation and drying times to suit the load. A very important feature for ensuring all bags and more complex instruments are thoroughly sterilised.

The large chamber size of 22 litre capacity enables more instruments to be sterilised at any one time. The Instaclave has space for 6 trays. 3 trays are included as well as your hose, drip tray and handle for ease of use. The door gasket can be replaced and the overheat cut out can be reset easily by the user, reducing your maintenance costs and down time.

### Why choose a non-vacuum Autoclave?

A non-vacuum autoclave is designed specifically for sterilisation of non-wrapped and solid instruments. It comes with the added benefit of being able to increase the sterilising and drying times to suit the load. There is also a larger version available.

### Integrated Safety Features

You will always feel comfortable whilst using the Instaclave due to its abundance of safety features. The manual and mechanical door locking mechanism ensures that it cannot be opened whilst sterilising due to the pressure forcing the second lock to engage. Overheating protection and a pressure safety valve are two other features that will guarantee you an easy and safe experience.

**Key Features:**

- Ability to increase sterilising and drying times
- 2 Cycle types
- Large chamber and tray size
- Safety features including manual and mechanical door locks and overheating protection
- Easy to use, clear buttons to switch between cycles, with clear lights to show your selection
- Non recycled water system

Technical Specifications	
Cycle Types	134°C Non Vacuum Cycle (drying cycle available) 121°C Non Vacuum Cycle (drying cycle available)
Chamber Size	22 litres
Tray Size	400mm x 180mm x18mm DWH
Capacity	6 Tray (Supplied With 3)
Dimensions	540mm x 500mm x 410mm WDH
Weight	45 Kg
Cycle Times	134°C Non Vacuum Cycle (30 or 40mins depending on drying cycle)* 121°C Non Vacuum Cycle (30 or 40mins depending on drying cycle)*
Power Requirements	220/240v, Current: 9 amps

**ADDITIONAL INFORMATION**

<b>SKU</b>	2545-V S3 A
<b>Overall Product Dimensions (mm)</b>	540 x 500 x 410 WDH
<b>Generic Size</b>	Medium
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Number of Cycles</b>	2
<b>Vacuum</b>	No
<b>Chamber Size</b>	22 Litres
<b>Autoclave Accessory/Spares</b>	Yes

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights