

BURTONS AIR BREAK ASSEMBLY - CLEARANCE

103-147-CL

**PRODUCT DESCRIPTION****Burtons Air Break Assembly**

Each anaesthetic machine will require an air break for use with the AGSS.

The air break eliminates any induced flow from the patient circuit by the negative pressure of the scavenging system.

Unlike a Barnsley receiver, our air break has a visual flow indicator, an important safety feature since the flow of air is a relatively quiet operation and the pump is mounted some distance from the point of use. Air breaks can be wall mounted in close proximity to the AGSS outlets if you choose to have wall mounted anaesthetic machines or mounted on the frame of a mobile anaesthetic machines.

Air breaks comply with BS6834:1987.

S/N: B59419

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	103-147-cl

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights