

COCOON FORCED AIR WARMING SYSTEM WITH 3 REUSABLE WARMING BLANKETS

CWS5000-PROMO24



PRODUCT DESCRIPTION

Designed to be quiet, yet effective, this forced air-warming system heats blankets quickly.

Offer Includes Set of 3 Reusable Warming Blankets (Small, Medium & Large)

The Cocoon warming system offers a great user experience thanks to its lightweight, compact, and sturdy design.

The user-friendly unit is equipped with two fan speeds, two flow rate levels, and four temperature settings, all controlled by simple and clear push buttons on the front. Despite its quiet operation, the DC motor provides the blower with all the power it needs.

Combined with Cocoon blankets, it provides quick, effective warming without blowing excess air around the surgical site. The blankets are made from soft, comfortable, and lightweight non-woven fabric, that is completely latex-free, and its unique perforations provide uniform air distribution without the risk of residue.

Universal inlet port ties on each blanket ensure that the hose remains firmly in place and the strong fabric allows the use of the same blanket on the patient, from pre-warming, through surgery, to recovery. The blankets are designed for single use and can be quickly disposed of after surgery. They are also easy to store and transport, making them the ideal solution for any environment. Some blankets have multiple inlet port locations for the best access. (See below for details on each blanket).

Key Features of the Cocoon warming controller.

- 2 fan speeds and 2 levels of flow rates
- 4 temperature options

https://burtonsveterinary.com/

- Uses a DC motor for a quiet, effective blower.
- Controlled by robust software.
- Unit is light, compact, and sturdy.
- Very simple, user-friendly design
- Alarm lights for system fault notifications
- Washable protective hose cover

Page 1 Sales: 01622 834300 | Service: 01622 834350



- Built-in safety monitoring system
- Retractable hooks at the back
- IV pole bracket attached.
- Extra-long power cable to get as close to your patient as possible.

Key Features of Cocoon warming blankets.

- Comfortable and lightweight fabric
- Microscopic perforations for quick warming without excess air being blown around the surgical site.
- Universal inlet port can be used with other blankets.
- Although these are single-use blankets, the fabric is strong enough to use multiple times throughout the patient's treatment, from surgery to recovery.
- Universal ties help to secure the blankets to the blower.

Usage Guide:

- An incontinence pad may be placed under the animal in the central area.
- Place the blanket either under or over depending on the blanket model.
- White side of the blanket faces the animal.
- Insert the heater hose into any blanket inlet port.
- Don't occlude the inflation tube at the blanket inlet or when closing a cage door.
- Monitor animal temperature frequently.
- Replace blanket if wet or soiled.
- Position the animal either dorsal or lateral.

ADDITIONAL INFORMATION

SKU	CWS5000-PROMO24
Animal	Mixed

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

https://burtonsveterinary.com/ Page 2 Sales: 01622 834300 | Service: 01622 834350