

MEDICAPTURE 170 VIDEO ARCHIVING UNIT

END-MCUSB170



Category: Image Capture & Management > Archiving Units > Flexible Endoscopy

Price: £2,590.88
Date Printed: 16/04/24

PRODUCT DESCRIPTION

The Medicap USB170 is the industry's easiest, high spec recording solution

Your images are saved directly onto USB flash drive or optional hard drive – a much simpler method than trying to capture images with a DVD recorder.

Key Features:

- **Medical Image Recorder** - The MediCap is even simpler than a normal camera. It records high quality still images from your endoscope, arthroscope, surgical camera, ultrasound etc.
- **Intuitive Use** - It is so easy to use, that no training is necessary. Just plug it in and go. You will be saving images within minutes of opening the box.
- **Medical grade** - With its steel case and fluid resistant front panel, the USB170 is built to withstand the rigors of the OR, and it is so compact it will fit on any cart.
- **USB hard drive** - perfect for archiving thousands of image on a single device
- **USB Keyboard** - Input patient data to organise your image library
- **Network Adapter** - Turn your USB170 into a mapped drive on your server
- **Foot pedal** - connect your USB170 to your scope's trigger or to a foot pedal

Specifications

- **Device Type:** Medical Image Recorder Records still images from standard medical video sources to a USB flash drive or optional USB hard drive
- **Case:** Rugged metal case with high-impact plastic faceplate
- **LCD Monitor:** 2 1/2", colour (for video confirmation only, not for diagnosis)
- **Buttons:** Sealed membrane, fluid resistant

- **Video In/Out:** PAL or NTSC format using s-video or composite BNC
- **USB Support:** USB2.0, high speed, 1 front connector, 3 rear connectors
- **Networking:** Ethernet T10/100 connection via Network Upgrade Kit (optional)
- **Footswitch/ trigger:** Hands-free capture Standard 3.5mm mini jack
- **Mechanical: Size:** 240 x 200 x 63mm / 9.5 x 7.9 x 2.5" Weight: 1.7kg / 3lb 12oz (appx.)
- **Environmental:** Storage: -40° to +85° C Operating: -20° to +40° C
- **Power:** AC 100-240V 50/60Hz 3-prong AC jack
- **Media:** USB flash drives or external USB hard drive
- **Image Formats:** JPEG, TIFF, PNG, DICOM
- **Image Size:** 1280 x 1024, 1024 x 768, 800 x 600, 640 x 480 pixels
- **Certifications:** CE 93/42/EEC, IEC 60601-1:2005, ANSI/AAMI ES60601-1:2005, CAN/CSA C22.2 No.60601-1:2008, E361228, FCC Part 15, FDA Class 1, RoHS, PSE, GOST-R
- **Kit contents:** MediCap USB170, USB flash drive, video cables, power cable, user guide

*The MediCap USB170 records still images only. If you are looking for a recorder that records still images and video, you may be interested in the USB200 (standard definition) or the USB300 (high definition and standard definition).

ADDITIONAL INFORMATION

SKU	END-MCUSB170
Animal	Mixed
Size	Small
Colour/Style	White

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights