

## OVAL LOOP RETRIEVAL SNARE 1.8MM X 180CM 2.5CM LOOP

END-F016



### PRODUCT DESCRIPTION

Reusable Oval Loop Grabber Snare With Inner Barbs 1.8mm x 180cm 2.5cm Loop (Hot Biopsy)

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	END-F016
<b>Channel Size</b>	1.8mm
<b>Autoclavable</b>	Yes
<b>Working Length</b>	1800mm

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

