

OVAL LOOP GRABBER SNARE 2.3MM X 240CM 3.0CM LOOP

END-F027

**PRODUCT DESCRIPTION****Oval loop grabber snare with inner barbs for optimum grip for foreign body retrieval**

- For use with 2.3mm channels or larger
- Ideal for foreign body retrieval in flexible endoscopes

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-F027
Channel Size	2.3mm
Working Length	2300mm

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights