

## BURTONS LAPAROSCOPIC POSITIONER

END-LAPP5



### PRODUCT DESCRIPTION

---

The new endoscopic positioner allows you complete control of your procedure with better visualization than ever before

Rotate your patient securely, left or right, without delays in surgical time and without disruption of the sterile field using the 3 angles that it offers. The rotation causes major abdominal organs to fall away from the body wall, allowing a better view. You are guaranteed to get the best position whilst keeping the patient safe and secure with multiple restraints and tie down points.

The anti-slip positioner fits securely to your table (fits variable table lengths) and with its robust, stainless steel frame, you are guaranteed an extremely high quality, long lasting addition to your practice.

Remove to clean.

#### Key Features:

- Facilitates quick rotation of patient without disrupting the sterile surgical field
- Provides easy visualization of the dorsal abdominal organs
- Versatile design attaches to existing surgery tables and accommodates any size patient
- Multiple tie-down and restraint points.
- Allows 3 different angles.
- High-quality stainless steel frame
- Easy to remove and clean.
- Adapts to various table lengths.
- Extra secure (no movement).
- Dual medical grade, high-quality rams.
- Includes patient tie down straps, various hooks, and spaces to adjust straps and ties, secure adjustments to fix positioner to the table (no need to purchase fixation clips).

**Only for use with a Burtons Table**

**Designed specifically to make these procedures faster and easier:**

Laparoscopy: Ovariohysterectomy (OHE), Ovariectomy, Cryptorchid, Abdominal exploration, Kidney biopsies.

Thoracoscopy: Pericardial sections, Pericardial windows, PD.

**Specifications:**

- Secure Index lock with spring pin return - Pull to release and place into new position slot.
- Frame hooks over ends of the table and are secured and adjusted using screw locks (can be seen in picture).
- Positioner will move 23 degrees from upright to flat.
- Table width: 600mm
- Fits standard and extended table tops: Minimum table top with standard clips 1182mm - Maximum table top with standard clips 1460mm.

**(Pro-Lift Table Not Included)**

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	END-LAPP5
<b>Brand</b>	Burtons

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights