

RIGID TELESCOPE FOR AI IN SMALL ANIMALS & RUMINANTS

END-SET-KB-41



Price: £4,590.92
Date Printed: 29/03/24

Category: Flexible Endoscopy > Telescopes > Endoscopes > Video Endoscopes

PRODUCT DESCRIPTION

AI Rigid Telescope & Instrument Set.

Key Features:

- 0 Degree Angle Of View
- 5mm Outer Diameter
- 310mm WL
- Autoclavable/ Soakable

Kit Includes:

- AI Rigid Telescope
- Transportation Case
- Trocars and Cannulas
- Atraumatic grasping Forceps
- Cleaning Brushes
- LED Light Handle

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	END-SET-KB-41
Telescope Degree (C)	0

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights