

SS-M SCOPESECURE MODULAR SYSTEM WALL MOUNTED SCOPE RACK

END-SS-M



PRODUCT DESCRIPTION

ScopeSecure Modular System Wall Mounted Scope Rack

- Holds one endoscope per rack
- Space saving design - purchase multiple racks to design a multi-scope storage system
- Made of durable, washable polycarbonate plastic
- Includes instructions and mounting hardware
- Width: 3 3/4 (3.75) in. / 9.5 cm
- Height: 4 1/4 (4.25) in. / 11.0 cm
- Depth (hook to base): 5 3/16 (5.1875) in. / 13.3 cm
- Material thickness: 1/5 (0.20) in. / 5.0 mm

ADDITIONAL INFORMATION

SKU

END-SS-M

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

