

FLEXIBLE SHORT OVAL CUP FORCEPS - 1.8MM

END-101-OCBUNDLE



Category: Instruments > Flexible Endoscopy

Price: £125.26
Date Printed: 29/03/24

PRODUCT DESCRIPTION

Flexible Forceps and Scissors - Short Oval Cup- Various Sizes Available

Highly reliable and easy to maintain, our flexible endoscopic instruments offer consistent performance and are CE and FDA approved.

- Constructed of special corrosion resistant steel alloys
- Autoclavable handles
- Colour coded and marked for identification

PRODUCT OPTIONS & ADD ONS

Working Length

1x Short Oval Cup 1.8mm D x 1750mm L (END-1015)
1x Short Oval Cup 1.8mm D x 1200mm L (END-1016)
1x Short Oval Cup 1.8mm D x 1900mm L (was £208.78)

ADDITIONAL INFORMATION

Animal

Mixed

SKU

end-101-OCbundle

Colour/Style	Stainless Steel
Catagory Sub Heading	Range of Working Lengths Available
Channel Size	1.8mm
Autoclaveable	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights