

STRAIGHT AIRWAY ADAPTOR (10PK)

V1100



Category: Patient Monitoring > Monitor Spares

Price: £173.25

Date Printed: 17/04/24

PRODUCT DESCRIPTION

Straight airway adaptor

- Airway Adaptor 13 ID x 15 OD mm.
- For monitors: V3395, V3404P, V8200, V8401,V8402, V90041, V9400, Advisor Monitor and most other side stream CO2 Monitors with a luer connector
- Sold in packs of 10

ADDITIONAL INFORMATION

Animal	N/A
SKU	V1100
Brand	SurgiVet
Type	Adaptor
Monitor Fit	Surgi Vet Advisor

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights