

## ECG MONITOR LEAD SETS

V3110TL-CONFIG\*



### PRODUCT DESCRIPTION

---

## ECG monitor lead sets with crocodile clips

- For monitors V3404, V3404P, Advisor
- Available in 3 lead 40" or 60"

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	v3110tl-config*
<b>Brand</b>	SurgiVet
<b>Type</b>	Lead
<b>Monitor Fit</b>	Surgi Vet Advisor

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights