

## MINDRAY ECG TRUNK CABLE

IMEC-0010-30-42719



### PRODUCT DESCRIPTION

---

## ECG Trunk Lead

- 3/5 Lead
- Adult or pediatric use
- 12 Pin
- Defibrillator proof
- AHA / IEC

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	IMEC-0010-30-42719
<b>Brand</b>	Mindray
<b>Type</b>	Lead, Other
<b>Monitor Fit</b>	Mindray iMEC 8, Mindray iPM12

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights