

END-MVE-9215 BRONCHSCOPE - BIOPSY CHANNEL: 2.8MM - DISTAL-END-9.2MM

END-MVE-9215



PRODUCT DESCRIPTION

The complete 9215 bronchoscope system has a broad range of applications and is the perfect solution for in surgery or ambulatory diagnosis

The MVE 9215 is a self-contained ambulatory endoscopy system. Replace your old fibrescope today with the latest, portable videoscope. With an integral screen, battery and recording capabilities, this bronchoscope is suitable for a broad range of applications such as; ear, nose, throat, bladder and urinary tract.

Show your images to clients and record images direct to SD card

For large animal and equine use.

Key Features:

- Working length 1500mm
- Diameter of distal-end: 9.2mm
- Diameter of instrument channel: 2.8mm
- Range of bending: U210 D90 L100 R100
- 3.5 inch LCD Monitor
- LED portable light source
- One case in car
- 8GB SD Card

Supplied With:

- Leak tester
- Biopsy forceps
- Battery pack
- Video recording head for moving image and stills
- Cleaning brush
- Aspiration catheter
- Operation manual

- Channel irrigator
- Protective case for transport
- One year free service and maintenance contract

ADDITIONAL INFORMATION

| | |
|-----------------------------------|------------------------------------|
| Animal | Equine |
| SKU | END-MVE-9215 |
| Colour/Style | Black |
| Brand | Aohua |
| Biopsy Channel | 2.8mm |
| Working Length (mm) | 1500 |
| Angulations | Up:210/ Down:90/Left:100/Right:100 |
| Fitting | Universal |
| Battery Operated/ Portable | Yes |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights