

## 8525- BURTONS MIXED PRACTICE VIDEOSCOPE

END-KW8525



### PRODUCT DESCRIPTION

---

## Burtons first fully immersible 8.0mm Videoscope - ideal for mixed practice

Ideal for small animal upper and lower Gi tract and some URT in medium to large breed dogs, aquariums, zoos, URT in equine, guttural pouches, and cystoscopy

#### Key Features:

- 8.5mm outer diameter
- 1500mm working length
- 2.5 mm biopsy channel
- 4-way angulations 210 degrees up / 100 degrees down
- 90 degree left/right

#### Supplied in transportation case with:

- Biopsy forceps x 2
- Cleaning brush x 2
- Flushing adapters
- Water bottle
- Spare piston kit
- Spare biopsy cap
- Waterproof cap

- 3 x 2 viewing and recording screen
- Video cables
- SD card
- User Manual

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	END-KW8525
<b>Colour/Style</b>	Black
<b>Brand</b>	Burtons
<b>Biopsy Channel</b>	2.5mm
<b>Working Length (mm)</b>	1500
<b>Angulations</b>	4 way
<b>Outer Diameter (mm)</b>	8.5

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights