

## TEMPOMATIC INFRA-RED TAP

115-267



### PRODUCT DESCRIPTION

---

## No contact with the tap means no risk of transmission!

With built in infrared spout sensor and automatic duty flush (for ~45 seconds 24 hours after the last use) this infra red tap is the perfect option for scrub sinks. As there is no contact made with the tap so hygiene is an absolute cinch!

The tap is easy to install due to it not requiring a separate control box or additional costs to provide an electrical power supply, unlike other infra-red taps. Instead this tap is powered by a long life lithium battery; typically lasting 3-5 years.

An automatic shut off and regulated flow at 6l/min (instead of the average 12~24l/min without auto shut off) means that you will save a substantial amount of water compared to conventional taps.

#### Features

- For temperature control this system can also be plumbed with a thermostatic mixer valve (stock code 115-268).
- Connection size: 1/2" BSP.(m)
- Supply pressure 0.5~5bar.
- Infra-red presence detection
- Flow restrainer set at 3 lpm at 3 bar
- M1/2" solenoid valve with filter
- Automatic duty flush (for ~45 seconds 24 hours after the last use)

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	115-267
<b>Colour/Style</b>	Stainless Steel

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights