


**SURGEONS LEVER WALL MOUNTED MIXER TAP**

115-261

**PRODUCT DESCRIPTION**

**Finished in a high shine chrome, this high quality mixer tap is designed for the riggers of commercial use**

**Features**

- Large 6" levers,
- ¼ turn ceramic disc for rapid on/off flow
- Swivelling spout and divided flow

**ADDITIONAL INFORMATION**

|                     |                 |
|---------------------|-----------------|
| <b>Animal</b>       | N/A             |
| <b>SKU</b>          | 115-261         |
| <b>Colour/Style</b> | Stainless Steel |
| <b>Operation</b>    | Surgeon (Lever) |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights