

**RIGID GOOSE 6" NECK SPOUT**

115-256-S

**PRODUCT DESCRIPTION**

**Finished in a high shine chrome, this high quality 6" gooseneck spout is designed for use with knee or foot pedal valves for a totally hands free start and stop**

**Key Features:**

- Adjustable gooseneck spout design
- Leak proof double O-ring spout design
- Stainless steel spout

**Specifications:**

- 6" Spout
- Rear mounted
- Requires separate knee (115-265) or foot (115-259) pedal valve
- Operating Temperature: 40-180°F (5-83°C)
- Operating Pressure: 15-125psi

**ADDITIONAL INFORMATION****Animal**

N/A

**SKU**

115-256-S

**Colour/Style**

Stainless Steel

**Operation**

Hand, Knee

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights