

## 20L WASTE BIN WHITE LID (FOR GENERAL USE)

225-458-W



### PRODUCT DESCRIPTION

---

Waste Bin White Lid (For General Use) 20L

Flame retardant construction, compliant with the new HTM07-01 (Safe Management of Healthcare Waste) standard set by the NHS.

Sizes HxWxD: 20 litre = 51x32x33cm

Product Features

- Hands free operation, foot pedal must be used to open
- Silent – dampened closing lid as standard
- Sealed base prevents leakages and facilitates easy cleaning
- Integrated wall mounting bracket as standard
- Corrosion resistant galvanised bodies with stainless steel foot pedal
- Designed and tested to HTM 83
- Elasticated cord hold bags firmly in place while the nylon loop makes changing bags easier

### ADDITIONAL INFORMATION

<b>SKU</b>	225-458-W
<b>Animal</b>	N/A

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights