

## FLAT OPERATING TABLE TOP

225-150-BUNDLE



## PRODUCT DESCRIPTION

---

### Standard flat table top

Standard flat top constructed of high quality stainless steel. Welded and polished corners allow the top to be thoroughly wiped clean with no hidden dirt traps. Complete with tie down points to secure the patient.

#### Why not make your top Heated ?

The heated top option available on a range of table top assists in maintaining body temperature during use. The temperature is controlled via an easy to use mounted panel, simply select the maximum temperature you want the table top to reach and the sensor within the table will then regulate the temperature once reached

#### Lead Lining

Available on a range of tops, lead lined tops are A high quality stainless steel top with built in lead lining, an Ideal combination for X-ray procedures and Operations/ Examinations where space is at a premium in the practice.

#### Don't forget to protect your top with a Burtons table mat!

We use a 2.5mm thick black matting made from a manmade material that is non-marking and has excellent chemical and abrasive resistance. Mats are provided with stainless steel securing clips for the two short ends, these help hold the mat in place on tables without a recess for the mat to lay in.

## Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

## Top Specifications

- Standard Top Size: 122 x 60 cm/ Heated Top Size: 137 x 60cm
- Heated Option
- Leadlined Option

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	225-150-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights