

## COMPOSITE RUBBER MAT WITH CLIPS FOR FLAT TOP TABLE

225-100-BUNDLE



### PRODUCT DESCRIPTION

## Protect your top with a Burtons table mat

We use a 2.5mm thick black matting made from a rubber composite material that is non-marking and has excellent chemical and abrasive resistance, perfect for use with high-level disinfectants. Mats are provided with stainless steel securing clips for the two short ends, these help hold the mat in place on tables without a recess for the mat to lay in.

- Standard top black mat With clips: 600 X 1220mm
- Extended top black mat with clips: 600 X 1370mm
- Custom Sizes available - Please contact our team.

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	225-100-bundle
<b>Colour/Style</b>	Black
<b>Brand</b>	Burtons

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights