

MGB ML-L ARISTO LIGHT SOURCE

END-865-00400A



Category: Flexible Endoscopy > Light Sources

Price: £6,414.54
Date Printed: 29/03/24

PRODUCT DESCRIPTION

MGB ML-L 300 Watt Xenon light source with integrated emergency halogen back up lamp - providing natural illumination

With L3 Laser/ LED hybrid technology with 7 inch touchscreen.

Key Features:

- Natural illumination
- Quick lamp exchange system
- Universal light guide connector
- 3 years manufactures warranty

Technical specification:

- Weight 14kg
- Dimensions 395 x 160 x 340 mm
- Bulb type 300-watt Xenon / 150-watt halogen emergency lamp
- Colour Temp approx. 6000 K
- Universal light guide connection for Xenon / AL3 type MGB/Wolf

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-865-00400A
Colour/Style	Stainless Steel
Brand	MGB
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Suitable For	MIS, Rigid Telescopes
Illumination Type	L.E.D
Fitting	Universal

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights