

## MAYO TROLLEY

115-06-BUNDLE



## PRODUCT DESCRIPTION

---

### Ideal for theatres where space is at a premium or simply to reduce clutter

With a sturdy Stainless steel construction, this mayo trolley with its removable tray and adjustable height is the perfect addition to your operating room. Lockable wheels ensure easy placement and removable tray allows for quick set up and cleaning.

Available with an optional bowl and 38in transfusion pole.

#### Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high-quality stainless steel materials this product has been engineered from the ground up to last!

#### Specifications

- Adjustable Tray Height: 88-130cm.
- Tray size: 51cm x 38cm
- 50mm locking castors

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	115-06-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights