

SERRATED CUP FORCEPS - 1.8MM

END-10-SF



Price: £0.00

Date Printed: 17/04/24

Category: Instruments > Flexible Endoscopy

PRODUCT DESCRIPTION

Flexible Forceps and Scissors - Serrated Cup - Various Sizes Available

Highly reliable and easy to maintain, our flexible endoscopic instruments offer consistent performance and are CE and FDA approved.

- Constructed of special corrosion resistant steel alloys
- Autoclavable handles

PRODUCT OPTIONS & ADD ONS

Serrated Cup Forceps

- 1x Serrated Biopsy Forceps 1.8mm D x 1200mm L (END-1024)
- 1x Serrated Biopsy Forceps 1.8mm D x 1750mm L (END-1064)
- 1x Serrated Biopsy Forceps 1.8mm D x 1900mm L (END-1094)

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-10-SF

Category Sub Heading**Procedure Type****Channel Size**

Range of Working Lengths Available

Laparoscopy, MIS

1.8mm

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights