

2 PRONG FORK GRASPERS - 1.8MM

END-10-2PF



Price: £0.00
Date Printed: 17/04/24

Category: Instruments > Flexible Endoscopy

PRODUCT DESCRIPTION

1.8mm 2 prong flexible fork graspers

PRODUCT OPTIONS & ADD ONS

Fork Sizes

- 1x 2 Prong Fork Grasper - 1.8mmD x 1200mmL (END-1026)
- 1x 2 Prong Fork Grasper1.8mm D x 1750mm L (END-1060)
- 1x 2 Prong Fork Grasper 1.8mm D x 1900mm L (END-1080)

ADDITIONAL INFORMATION

Animal	N/A
SKU	END-10-2PF
Catagory Sub Heading	Range of Working Lengths Available
Procedure Type	Bronchoscopy
Channel Size	1.8mm
Autoclaveable	Yes

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights