

SINGLE USE AND DISPOSABLE BIOPSY FORCEPS 1.8MM - 2.3MM

END-DISPOSABLE-FORCEPS-2.3



PRODUCT DESCRIPTION

A range of sterile single use and disposable biopsy forceps for scopes with 1.8mm and 2.3mm channels, suitable for small animal and equine

ADDITIONAL INFORMATION

| | |
|-----------------------------|------------------------------------|
| Animal | Mixed |
| SKU | end-disposable-forceps-2.3 |
| Catagory Sub Heading | Range of Working Lengths Available |
| Channel Size | 2.3mm |
| Autoclavable | No |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

