

DISPOSABLE BIOPSY FORCEPS - 2.3MM

END-DISPOSABLE-FORCEPS-2.3



Price: £0.00

Date Printed: 29/03/24

Category: Instruments > Flexible Endoscopy > Diagnostic Imaging

PRODUCT DESCRIPTION

2.3mm disposable biopsy forceps for small animal and equine scopes with 2.3mm channels and larger

PRODUCT OPTIONS & ADD ONS

Choose Your Length

- 1x Box of 10 coated disposable oval cup forceps 2.3mm x 2300mm
- 1x Box of 5 Disposable Oval Cup Forceps 2.3 x 2100mm (END-OJ-2421D)

Equine Specific

- 1x Box of 5 Coated Disposable Oval Cup Forceps 2.3 x 3500mm (END-OJ-2435D)

ADDITIONAL INFORMATION

Animal	Mixed
SKU	end-disposable-forceps-2.3
Catagory Sub Heading	Range of Working Lengths Available
Procedure Type	Bronchoscopy
Channel Size	2.3mm

Autoclaveable

No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights