

CLEANING BRUSHES - 2.3MM

END-2.3-CLEANING



Category: Instruments > Flexible Endoscopy

Price: £0.00
Date Printed: 16/04/24

PRODUCT DESCRIPTION

Variety of cleaning brushes for all 2.3mm channels

PRODUCT OPTIONS & ADD ONS

Small Animal

1x Autoclavable cleaning brush 2.3mm x 2100mm (END-CY-2421D)
1x Autoclavable Cleaning Brush 2.3mm x 3500mm (END-CY-2435D)

Equine

1x Autoclavable Cleaning Brushes 2.3mm X 4000mm (END-CY-2440A)

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-2.3-CLEANING
Catagory Sub Heading	Range of Working Lengths Available
Procedure Type	Arthroscopy
Channel Size	2.3mm

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights