

## SERRATED CUP FORCEPS - 2.3MM

END-SERRATED-FORCEPS-2.3



Category: Instruments > Flexible Endoscopy

Price: £0.00

Date Printed: 17/04/24

### PRODUCT DESCRIPTION

**Serrated cup forceps designed for small animal and equine scopes with 2.3mm channels or larger.**

- Autoclavable
- Various working lengths available

### PRODUCT OPTIONS & ADD ONS

#### Choose your length

- 1x Serrated Biopsy Forceps 2.3mm D x 1200mm L (END-1234)
- 1x Serrated Biopsy Forceps 2.3mm D x 1650mm L (END-1264)
- 1x Serrated Biopsy Forceps 2.3mm D x 2350mm L (END-1294)

#### Equine Specific

- 1x Serrated Biopsy Forceps 2.3mm D x 3500mm L (END-3794)

### ADDITIONAL INFORMATION

**Animal**

Mixed

**SKU**

END-SERRATED-FORCEPS-2.3

**Category Sub Heading****Procedure Type****Channel Size****Autoclaveable****Working Length**

Range of Working Lengths Available

Bronchoscopy

2.3mm

Yes

1200mm

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights