

**SERRATED CUP FORCEPS WITH SPIKE -2.3MM**

END-SERRATED-2.3-SPIKE

**PRODUCT DESCRIPTION****Serrated Biopsy Forceps with a spike for 2.3mm endoscopy channels and larger**

- Autoclavable
- Various working lengths available
- With spike

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	END-SERRATED-2.3-SPIKE
<b>Catagory Sub Heading</b>	Range of Working Lengths Available
<b>Channel Size</b>	2.3mm
<b>Autoclavable</b>	Yes

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights