

### 3 NAIL GRASPERS - 2.3MM

END-NAIL-GRASPER-2.3



## PRODUCT DESCRIPTION

### 2.3mm Nail Grasper for foreign body retrieval in flexible endoscopy procedures

- Autoclavable
- Various lengths available
- For 2.3mm channels or larger

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	END-NAIL-GRASPER-2.3
<b>Catagory Sub Heading</b>	Range of Working Lengths Available
<b>Channel Size</b>	2.3mm
<b>Autoclaveable</b>	Yes

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

