

OVAL SNARES - 2.3MM

END-OVAL-SNARE-2.3



Price: £0.00
Date Printed: 28/03/24

Category: Instruments > Flexible Endoscopy

PRODUCT DESCRIPTION

Oval Snares for foreign body retrieval in various sizes

For use with 2.3mm channels and larger

PRODUCT OPTIONS & ADD ONS

Loop Size

- 1x 2.3mm Oval Snare - 30mm x 2350mm L (END-1100)
- 1x 2.3mm Oval Snare 25mm - 2350mm L (END-1101)
- 1x 2.3mm Oval Snare Micro - 20mm x 2350mm L (END-1102)

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-OVAL-SNARE-2.3
Catagory Sub Heading	Range of Working Lengths Available
Procedure Type	Gastroscopy
Channel Size	2.3mm
Working Length	2350mm

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights