

NSK HIGH SPEED HANDPIECE

148-550



PRODUCT DESCRIPTION

High-speed handpiece with a patented clean head system, designed to prevent oral fluid or foreign particles from entering into the handpiece head.

The system is activated instantaneously in response to air pressure change as the handpiece stops spinning.

Key Features:

- Ceramic Bearings - long lasting
- Clean Head System for enhanced infection control
- Ultra Push Chuck - simple yet secure fastening
- Single Spray
- Speed : 380,000~480,000 min-1
- Head Size: ø10.3 x H13.8 mm
- Midwest 4 hole
- Single Spray - efficient cooling
- Ergo Grip - secure grip for your confidence
- Handpieces are autoclavable to 135°.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	148-550
Generic Size	X-Small
Examination Type	Oral Examination

Fitting Type

4 Hole

Burr Fit

F/G

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights