

## CHEEK DILATOR - EXTRA SMALL

128-006-W



### PRODUCT DESCRIPTION

## Designed for small rabbits, chinchillas, and other small rodents

Usually under approx. 1.5kg body weight. A cheek dilator retracts cheeks to allow visualization and access to teeth.

#### Key Features:

- Lightweight and takes up less space than solid paddle design
- Spring loaded for quick and safe use
- Stainless steel
- Use alongside mouth gag
- Wire frame paddles fit snugly into cheek surface - take up less space compared to a solid paddle type

### ADDITIONAL INFORMATION

|                               |              |
|-------------------------------|--------------|
| <b>Animal</b>                 | Small Furies |
| <b>SKU</b>                    | 128-006-W    |
| <b>Generic Size</b>           | XX-Small     |
| <b>Dental Instrument Type</b> | Dilators     |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights