

**STRAIGHT HANDPIECE**

148-400

**PRODUCT DESCRIPTION****Straight Handpiece for use with HP style burs, ideal for rabbit dentals****Key Features:**

- For use with HP style burs
- Ideal for rabbit dental for 'E' fitting handpieces
- Suitable for use with micro motor systems and slow speed air motors
- Handpieces are autoclavable to 135°C
- Handpieces come with 6 months warranty

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	148-400
<b>Examination Type</b>	Oral Examination, Small Furries
<b>Fitting Type</b>	E
<b>Burr Fit</b>	HP Shank

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights