

STRAIGHT HANDPIECE

148-400



PRODUCT DESCRIPTION

Straight Handpiece for use with HP style burs, ideal for rabbit dentals

Key Features:

- For use with HP style burs
- Ideal for rabbit dental for 'E' fitting handpieces
- Suitable for use with micro motor systems and slow speed air motors
- Handpieces are autoclavable to 135°C
- Handpieces come with 6 months warranty

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	148-400
Autoclaveable	Yes
Examination Type	Oral Examination, Small Furrries
Fitting Type	E
Burr Fit	HP Shank

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights