


**EMS SCALER TIP PS**

118-218

**PRODUCT DESCRIPTION****Removes 95% of all calculus & mineralized deposits****Key Features:**

- High-quality surgical-grade stainless steel
- Intricate manufacturing process with up to 31 production steps
- Jewel-grade polishing
- Absolutely linear oscillations for unequalled precision and predictable treatment outcome
- Bactericidal action through cavitation and acoustic streaming
- Gentle treatment for maximum conservation of tooth substance
- Unique combination of torque wrench and instrument holder, significantly reducing the risk of instrument loss and contamination

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	118-218
<b>Brand</b>	EMS
<b>Autoclaveable</b>	Yes

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights