

## URINARY TIPS

118-2-URINARYTIPS



## PRODUCT DESCRIPTION

**All tips are supplied with the CombiTorque mounting system which is used for proper tip mounting and secure storage of each individual instrument.**

Your Piezon S scaler can be fitted with these uniquely developed urinary tips, which are inserted into the urethra to pulverize wedged calculus, a common problem in cats.

- Dog Tip Size:11.4mm
- Scaler tips are autoclavable to 135°.

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	118-2-urinarytips
<b>Category Sub Heading</b>	Range of Sizes Available
<b>Autoclaveable</b>	Yes

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights