

## VISTA 40 STORAGE TROLLEY

225-420-BUNDLE40



### PRODUCT DESCRIPTION

---

## The Vista 40 is a compact medical/storage trolley ideal for under worktops

The range of Vista trolleys is one of the largest and most flexible on the market today, widely used by animal and human care professionals.

The Sunflower range of Vista medical trolleys and storage modules is extensive, with options to suit all applications and environments, from GP or hospital consultation and treatment rooms, to hospital wards and intensive care units. The product range is designed to be easy to clean, with smooth lines and minimal dirt traps. Each medical trolley or storage module can be accessorised with medi-rails, sharps box holders, glove box holders, alcohol gel holders, oxygen cylinder holders, laptop holders, and more.

#### Features

- Tilt & locking tray runner system
- Trays integrate with all other Vista products
- Low level (80cm), ideal for under worktops
- 100mm castors with front brakes
- Optional tamper resistant locking door
- Removable top shelf with integrated upstand
- Stainless steel handle
- Easy clean design and materials
- Supplied fully assembled
- 2 year warranty

**Dimensions:**

- Vista 40 Trolley 80cm x 56cm x 52cm HWD

**Trays:**

- Narrow/single depths 8cm x 31cm x 42cm HWD
- Narrow/double depths 16cm x 31cm 42cm HWD
- Wide/single depths 8cm x 47cm x 42cm HWD
- Wide/double depths 16cm x 47cm x 42cm HWD

**For accessories, see below. Please note, the following can only be fitted when ordered with the trolley:**

- Medirail
- Protective buffers
- Folding shelf
- Sharps box holder

**ADDITIONAL INFORMATION**

|  |                        |
|--|------------------------|
| <b>Animal</b>                          | Mixed                  |
| <b>SKU</b>                             | 225-420-bundle40       |
| <b>Overall Product Dimensions (mm)</b> | 80cm x 56cm x 52cm HWD |
| <b>Colour/Style</b>                    | White                  |
| <b>Brand</b>                           | Sunflower Medical      |
| <b>Kick Bowl</b>                       | No                     |

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights