

**QS12 - ULTRAWAVE CLEANING BATH - 12.5L**

116-619

**PRODUCT DESCRIPTION****Ultrawave's Q-series baths allow accurate cycle parameter control and effective cleaning performance for all applications**

Offering repeatable and traceable cleaning performance for a wide range of industries. The memory card port allows traceability of cycle number, time and temperature to be stored on your PC. Heated, and with easily operated digital time and temperature control, the QS range has tank capacity ranging from 4.5 to 25 litres and is supplied complete with a stainless steel basket and ABS plastic lid.

Please note: 3 year warranty included when you order with us\*

*\*1 year as standard, 2 years extra warranty applicable if registered online with Ultrawave within 30 days of placing your order*

|                          |  |
|--------------------------|--|
| Overall Dimensions       | 345 x 285 x 365 (W x D x H)  |
| Internal Tank Dimensions | 295 x 245 x 200mm  |
| Basket Dimensions        | 265 x 215 x 140mm  |
| Capacity                 | 12.5 Litres  |
| Constrution              | 304 Stainless Steel Tank   |
| Control                  | Digital Control<br>Time (0 - 99minutes)<br>Temperature (Ambient to 70°C) |
| Operating Feaxtures      | Memory Card Port For Digital Traceability<br>Water Level Sensor          |

|                     |                    |
|---------------------|--------------------|
| Ultrasonic Power    | 200W               |
| Operating Frequency | 32 to 38kHz        |
| Heating Power       | 300W               |
| Required Services   | 13amp 3 Pin Socket |

## ADDITIONAL INFORMATION

|  |                             |
|--|-----------------------------|
| <b>Animal</b>                          | Mixed                       |
| <b>SKU</b>                             | 116-619                     |
| <b>Overall Product Dimensions (mm)</b> | 345 x 285 x 365 (W x D x H) |
| <b>Brand</b>                           | Ultrawave                   |
| <b>Ultrasonic Bath Controls</b>        | Digital                     |
| <b>Ultrasonic Bath Litres</b>          | 12.5                        |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights