

QS18 - ULTRAWAVE CLEANING BATH - 17.5L

116-620



Category: Infection Control > Ultrasonic Baths > Ultrasonic Cleaners > Medical > Ultrasonic Baths

Price: £2,191.00
Date Printed: 04/10/23

PRODUCT DESCRIPTION

Ultrawave's Q-series baths allow accurate cycle parameter control and effective cleaning performance for all applications

Offering repeatable and traceable cleaning performance for a wide range of industries. The memory card port allows traceability of cycle number, time and temperature to be stored on your PC. Heated, and with easily operated digital time and temperature control, the QS range has tank capacity ranging from 4.5 to 25 litres and is supplied complete with a stainless steel basket and ABS plastic lid.

Overall Dimensions	345 x 285 x 365 (W x D x H)
Internal Tank Dimensions	330 x 300 x 200mm
Basket Dimensions	285 x 260 x 140mm
Capacity	17.5 Litres
Construction	304 Stainless Steel Tank
Control	Digital Control Time (0 - 99minutes) Temperature (Ambient to 70°C)
Operating Features	Memory Card Port For Digital Traceability Water Level Sensor
Ultrasonic Power	300W
Operating Frequency	32 to 38kHz

Heating Power	450W
Required Services	13amp 3 Pin Socket

ADDITIONAL INFORMATION

Animal	Mixed
SKU	116-620
Overall Product Dimensions (mm)	375 x 345 x 365 (W x D x H)
Brand	Ultrawave
Ultrasonic Bath Controls	Digital
Ultrasonic Bath Litres	17.5

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights