

ULTRACLEAN M2 X 6

116-601



Category: Infection Control > Consumables > Ultrasonic Baths > Ultrasonic Cleaners > Disinfectant > Medical > Ultrasonic Consumables

Price: £104.95
Date Printed: 16/04/24

PRODUCT DESCRIPTION

Ultraclean M2 is a general purpose cleaning detergent, suitable for use in a wide variety of applications.

Pack of 6

Ultraclean M2 can be used to remove contaminants from plastics, glass, metals and other materials.

Key Features:

- General purpose cleaning detergent
- Suitable for cleaning medical instruments
- Reduces surface tension
- Non-foaming

(If processing aluminum or other soft metals, Ultraclean M2 should not be used)

ADDITIONAL INFORMATION

Animal	Mixed
SKU	116-601
Brand	Ultrawave

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights