

RO WATER PURIFICATION SYSTEM

139-000



PRODUCT DESCRIPTION

Thanks to the deionization stage, this RO system will remove almost 100% of tap water contaminants

Key Features:

- Ideal water for your autoclave
- Reverse osmosis water on demand
- No storage of bulky containers
- Mount it anywhere
- Plumb into your existing mains water supply
- Simple to install in less than 1 hour
- Removes up to 99.8% of all dissolved liquids

The Osmotics RO system includes:

- 1 x 75 US gallon per day membrane
- 1 x 5 micron Sediment Pre-filter
- 1 x Carbon Block Pre-filter
- 1 x 700ml Deionisation Post-filter
- 1 x Feed Water Clamp
- 1 x Waste Water Clamp
- 1 x Filter Housing Wrench

- 2 Metres of Red Tubing
- 2 Metres of Green Tubing
- 2 Metres of Blue Tubing
- 1 x Membrane Flush Kit *

* Using the flush valve regularly increases membrane life.

Don't forget to buy a TDS meter which monitors the PPM (particales per million) content of the mains supply, filtered R/O water. These items can be found on the related products tab.

The effectiveness and productivity of RO systems increases with greater water pressure and temperature. The actual flow rate (in GPD) of any reverse osmosis system depends on your water pressure and temperature.

While manufacturers state that 60 psi is required at the membrane for full production, we know that not everyone's cold water is at 60 psi and 77 degrees F therefore, we remind you that you might get less than 75 GPD, perhaps 50% or even less if your pressure and temperature are very low.

Reverse osmosis self regulating booster pumps are available on the additional items below.

ADDITIONAL INFORMATION

SKU

139-000

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights