

## LABCOLD PHARMACY UNDER-COUNTER REFRIGERATOR 150 LITRE

125-041



### PRODUCT DESCRIPTION

**Powered by the unique Intellicold controller, this professional 150 litre pharmacy refrigerator has been designed from the ground up for the safe storage of cold chain medicines**

With superior temperature uniformity, low energy consumption and extra functions, this refrigerator is ideal for storing your sensitive pharmaceuticals in perfect conditions, helping busy professionals monitor regulations and requirements.

Labcold IntelliCold Pharmacy Refrigerator now with NEW touchscreen controller. Designed from the ground up for the safe storage of cold chain medicines as advised by Public Health England.

Featuring a 1 metre secondary temperature probe that can be placed anywhere inside the pharmacy fridge and a touchscreen controller which is completely flat so it can be thoroughly cleaned and disinfected. The controller displays alarms and records minimum and maximum temperatures for both temperature probes. It also records alarm states, fridge operation, door openings and temperature data for both probes at minute intervals.

You can only view and store this data on a computer by purchasing a specially encoded micro SD card. This data helps you with record keeping and regulatory inspections. Plus, if you have any issues, you can establish the exact duration and severity. You can find these in the related products below

#### Key Features

- Fitted with the NEW IntelliCold® touch screen controller
- 24/7 automatic temperature data logging which is downloadable via optional micro SD card

- Controller and alarm battery back up
- Door open alarm
- Power failure alarm
- High temperature alarm
- Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Access port
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Self-closing door
- Brackets for wall mounting
- Real time temperature display for two temperature probes
- Secondary (wandering) temperature probe
- Internal light
- Remote alarm contacts
- Optional micro SD card for transferring temperature data to your computer
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines

**Specifications:**

- Configuration Undercounter
- Capacity (Litres) 150 litres
- Height x Width x Depth 835mm x 600mm x 600mm
- Weight 51kg
- Shelves 3
- Temperature Range +2°C to +8°C
- Door Type Solid
- Lock Type Key Lock
- Temperature Display
- High/Low Temperature Alarm
- Door Open Alarm
- Mains Failure Alarm
- Max/Min Temp. Recording
- Alarm battery backup
- Auto Defrost
- Access Port
- Remote Alarm Contacts
- Integral Temp. data logging

**ADDITIONAL INFORMATION**

<b>SKU</b>	125-041
<b>Overall Product Dimensions (mm)</b>	835mm x 600mm x 600mm
<b>Colour/Style</b>	White
<b>Brand</b>	Labcold
<b>Alarm</b>	Yes
<b>Fridge Capacity</b>	150
<b>Door Type</b>	Regular

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights