

NEEDLE ELECTRODES

NEEDLE ELECTRODES



Price: £0.00

Date Printed: 28/03/24

Category: Electrosurgery > Accessories

PRODUCT DESCRIPTION

Needle electrodes: ('A' Series) - used mainly for incisions and excisions in precise surgical procedures

Also included in this section:

- Broad needle electrode, mainly used for fulguration and desiccation
- Vari-tip electrode, the tip can be pushed into shorted tip length for eyelid or skin incisions, the wire will eventually become brittle and can be extended and cut for continued use

PRODUCT OPTIONS & ADD ONS

Needle Electrodes

- 1x Broad Needle Electrode - Twin Pack
- 1x TA1 - Needle Electrode Ellman FFPF - Twin Pack
- 1x TA2 - Needle Electrode Ellman FFPF - Twin Pack
- 1x TA3 Needle Electrode Ellman FFPF - Twin Pack
- 1x Vari-Tip A8 With Spare Wire
- 1x Needle Electrode 3/8in - Twin Pack

ADDITIONAL INFORMATION

SKU

Needle Electrodes

Brand

Ellman

Animal

Mixed

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights