

TRIANGLE ELECTRODES

TC-TRIANGLE-ELECTRODES



Price: £0.00

Date Printed: 29/03/24

Category: Electrosurgery > Accessories

PRODUCT DESCRIPTION

Triangle electrodes: ('C' Series) - Used mainly for shaving lesions and contouring tissue.

PRODUCT OPTIONS & ADD ONS

Triangle Electrodes

- 1x C10 Triangle Electrode - Twin Pack
- 1x C11 Triangle Electrode
- 1x C12 Triangle Electrode - Twin Pack
- 1x C13 Triangle Electrode - Twin Pack
- 1x C14 Triangle Electrode - Twin Pack
- 1x C15 Triangle Electrode - Twin Pack
- 1x C7 Triangle Electrode - Twin Pack
- 1x C8 Triangle Electrode - Twin Pack
- 1x C9 Triangle Electrode - Twin Pack

ADDITIONAL INFORMATION

Animal	Mixed
SKU	TC-Triangle-electrodes
Brand	Ellman

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights