

## CONTROLLED DRUGS CABINET - 500W X 300D X 550HMM

115-493



### PRODUCT DESCRIPTION

---

## Single point locking controlled drugs cabinets that exceed the statutory requirement for the storage of controlled drugs

Ensuring security with a strong 1.5mm (16 gauge) mild steel, seamless welded construction and fitted with a BS3621 anti-pick & anti-drill, six-pin cylinder key retaining lock. Supplied with 2 working keys, two shelves and identity card.

#### Key Features:

- Supplied right hand hinging as standard, left hand available on request
- Single-point lock
- supplied with 2 Keys & ID Card
- 2 x Shelves
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs per profile
- Lighting transformer casing adds 105mm to the overall height
- 1.5mm (16 gauge) mild steel, seamless welded construction
- Mounting holes in back (fixings not included)
- Adjustable & removable shelves
- Modular design allows for easy interchange with HTM 63 storage systems
- Painted white finish
- Compatible with ABLOY® CLIQ™ Remote, integrating electronics into high-security mechanical disc cylinders & keys, combining electronic security & remote management in real time. Can be supplied with a micro-switch to allow connection to remote lighting/nurse call system by others (CD905), or supplied with a 12v DC LED lighting/warning system with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC

**Dimensions:**

- 500W x 300D x 550Hmm
- 20.0kg

**ADDITIONAL INFORMATION**

|  |                      |
|--|----------------------|
| <b>Animal</b>                          | N/A                  |
| <b>SKU</b>                             | 115-493              |
| <b>Overall Product Dimensions (mm)</b> | 500W x 300D x 550Hmm |
| <b>Colour/Style</b>                    | White                |
| <b>Brand</b>                           | Bristol Maid         |

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights