

BURTONS HIBISCRUB DISPENSER BRACKET

115-292-CONFIG

**PRODUCT DESCRIPTION**

High quality, strong, steel Hibiscrub dispenser bracket for mounting on walls or Burtons scrubs sinks (Additional bracket required)

- Large lever for easy dispensing
- Easily wall mounted
- Will not bend or snap with use
- Perfect match for Burtons scrub sinks with addition bracket (115-294)

*Scrub Sink Bracket and bottle not included

ADDITIONAL INFORMATION

Animal	N/A
SKU	115-292-config
Colour/Style	Silver
Brand	Burtons

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights