

SURGEONS DECK MOUNTED 6" LEVER TAP

115-250-3



PRODUCT DESCRIPTION

Ideal for use on surgical sinks

Finished in a high shine chrome, this high quality tap is fitted with 6" levers and ceramic ¼ turn seats for rapid on/off flow

Features

- Swivelling spout and divided flow
- Offset mountings enable fitting on any sink / surface with holes centres between 160mm and 240mm.
- Inlet size 1/2" BSP
- Deck Mounted

ADDITIONAL INFORMATION

| | |
|---------------------|-----------------|
| Animal | N/A |
| SKU | 115-250-3 |
| Colour/Style | Stainless Steel |
| Operation | Surgeon (Lever) |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights